Automatic Jaw Motion Tracking & Analysis Software User's Manual

2019.8 Rise Corporation

 $\mathbf{2}$

Contents

4
4
.11
.10
.12
.15
.15
.19
.21
.23
.24
.24
.26

1. Create Analysis Data

The software can analyze frontal, sagittal, and horizontal movements from Movie. Combining frontal and sagittal planes of analysis data can analyze horizontal and 3D data. However, in the case of compositing, it is necessary to match the simultaneous shooting and the analysis start frame.

In the case of 3D data, it can be synthesized with the average data obtained from frontal and sagittal analysis.

1-1. Analysis from Movie

Create patient analysis data.

Import the movie you shot in advance to WinCeph. Create analysis data of frontal plane.



In the image management screen, select the captured frontal movie and click [Jaw Motion].



Then select [Frontal] and click [OK]. (To create a [Sagittal] plane or [Horizontal] plane, just select each item on the left screen. The procedure is the same.) The [Marker Setting] dialog is displayed, displaying the first frame of the movie. Designate a marker that displays the trajectories of the upper and lower jaws.



Use the mouse drag, enclose and set the markers in order of upper and lower jaw. (The center of the marker is reflected in the trajectories of upper and lower jaws)

The upper right of the screen [Display Size], the larger the number, the screen will be larger, and it will be easier to specify the marker.

When setting is complete for both upper and lower jaws, click [OK].

- *When setting the marker, if it is difficult to specify the marker for the first frame image (due to blurring of the image, etc.), change the frame with the scroll bars at the bottom of the screen and specify the marker with clear frame of the image.
- *left rotation and right rotation are used; image analysis data can be created even when the camera is used vertically.
- *Mirror...Perform left-right reversal.

Next, the "Marker Setting / Tracking" dialog is displayed. Click [Start Tracking] to track the set marker. (For a detailed explanation, please refer to {2. Refer to Patient Analysis Data}) Start tracking and acquire markers for all upper and lower jaws for each movie.



When tracking is complete, the list on the left of the screen shows the frame number and the position of the upper and lower jaw markers on the screen.

Depending on the shooting conditions, markers may not be detected properly.

In this case, try switching the detection mode or changing the tolerance value.

When setting the detection mode to [Mode B], it is necessary to set the marker detection range.

* Setting of detection range



Specify the movement range of the upper and lower jaw markers by mouse drag with the [Set Detection Range].

Scroll bar and [<] [>] can change the frame of the movie.

When complete, click [OK].

After changing the detection mode or range setting, click [Tracking Start] again. When tracking is complete, click [NEXT Error] to check for marker acquisition errors.

Frames with errors (frames for which markers cannot be set) are displayed in red in the left list. (Position display is -1)

If there are many errors, change the detection settings etc. and try again.

In addition, if you want to analyze only a specific part of the movie, please specify [Start Frame] and [End Frame].



The error frame is displayed, and the marker that failed to be detected is displayed on the left of the screen.

Move the upper and lower jaw markers to the respective marker positions with mouse drag. Click [NEXT Error] again after manual marker setting. Repeat the manual marker setting until there are no errors. Click [OK] when there are no errors.



Calibration setting screen is displayed.



Make [Specify 2 Points] pressed. (It is pressed in the initial state) The mouse pointer changes to a green cross on the image, so drag the mouse according to the object whose length you can see (scale, marker, etc.).

You can also change the [Display Size] on this screen.

As it gets bigger the [View Size], you can more accurately "distance designation".

(If [Move Image] is pressed, the image can be moved by dragging the mouse.

It is used to scroll the screen when the image size is bigger.)

Click [OK] when [Specify 2 Points] is complete.

Designation of Distance Between Two Points $\qquad imes$	
Specify the Distance Between Two Points in mm.	Enter the length between the two specified points
0K Cancel	in "mm" and click [OK].

As the "Frontal" dialog is displayed, click [Save Data] to save the data.



The same analysis applies to sagittal planes and horizontal planes.

1-2. Synthesize Analysis Data

(1) Horizontal Composition

Using analyzed frontal and sagittal plane data, it is effective only when the animation of each plane is simultaneously shot and the analysis start is together.

. L Onen	🛃 Lat.	PA 🛞	Model 🛞	law Movement	
Open Analysis Window	Date ^	PhotoNo.	Age	C. Comment	
Show Info. Show Analytical Info.					
Ricketts View					
Copy Copy Analytical Info.					
∿r Jaw Movement Analysis					

Click [Jaw Motion] from the analysis list.

Select Process X
O Create 3D Data
Create Horizontal Data
🔿 Compare Analysis Data
OK Cancel

Select [Create Horizontal Data] and Click [OK]

A screen for selecting the data to be synthesized is displayed.

lorizontal Analysis			
** Use Only for Sim	ultaneous Sh	ooting	
Frontal Data			
Date	Age	Comment	Numbr
2019/7/30	39u2m		12
2003/ 7/29	23y2m		11
_	-		
<			>
Sagittal Data			
Date	Age	Comment	Numbe
2003/ 7/29	23y2m		11
<			>
		OK	Cancel

Check the frontal and sagittal data taken simultaneously, and click the [OK].



Analysis results of the combined horizontal plane are displayed.

(2) 3D View

Select Process	X Data Data	 Click [Jaw Motion] from the analysis list, select [Create 3D Data], and click [OK]. 								
OK Car	ncel									
Select 3D Constru	iction Da	ta			×					
Select 3D Constru	uction Met	hod								
Synthesize	3D Data (Only on Aver	age Cycle							
🔿 Combine 3D) Data wit	h Analysis Ef	fective Data	Specify	composit	ion type				
** Use Only Frontal Dat	for Simult	aneous Shoo	oting							
Date	-	Age	Comment		Numbe					
2019/	7/31	39y2m			11					
2019/	7/31	39y2m			12					
<					>					
Sagittal Da	ita									
Date		Age	Comment		Numbe					
2019/	7/31	39y2m			11					
<					>					
			OK	Can	cel					

- * "Synthesize 3D Data Only on Average Cycle" allows 3D compositing even if the frontal and sagittal movies are not shot simultaneously.
- * "Combine 3D Data with Analysis Effective Data" allows 3D compositing only if the frontal and sagittal videos are taken simultaneously and the analysis starts together.

(If the shooting is not simultaneous or the analysis start frame is different, the display of 3D composition is not performed correctly.)

Check any frontal and sagittal planes and click [OK].



* When combining only the average cycle

* In the case of 3D synthesis with analysis effective data



2. Refer to Patient Analysis Data

2-1. Refer to Analysis Data

Basic Information Images Mana	t agement Photo Gall	ery Anal	ysis List Cas	e Information	Slideshow	Select [Jaw Movement]
Open Open Analysis Window	Lat. Date ^ 2003/7/29 2019/7/30	PA OF Type Frontal Frontal	Model/ () PhotoN	aw Movement Age 23y 2m 39y 2m	Comment	
Show Analytical Info.	2003/ 7/29 2003/ 7/29 2003/ 7/29	Sagittal Horizontal 3D	0000000002 0000000000 0000000000	23y 2m 23y 2m 23y 2m		
Ricketts View Superimposition Ricketts						
Copy Copy Analytical Info.						
) Delete						

Select the data you want to refer from the analysis list.



(Click [Delete] to delete the selected data.)



Saves the current data. (marker setting, tracking setting, calibration, selection cycle, R rotation, L rotation, Mirror).

CalibrationCalibration

Enter the length between two specified points in mm unit and correct the distance.



Display the dialog of Cycle List.



Display the dialog of Speed Cycle.



..... Tracking

Make settings for markers and tracking again.



..... Display the entire cycle.

Marker setting and Tracking setting



Tracking Test...... Tests whether the marker position of the displayed image can be detected. The detected position is not retained.

Start Tracking...... Detects markers of all frames in the specified range.

Stop Stop marker detection.

Detection Mode

Mode A ... Normal mode

(Detected based on RGB values)

Mode B...Specified range detection mode

(Detected within the range specified by Set Deletion Range)

Mode C...Hue detection mode

(Detected based on the values of Hue, Chroma and Lightness)



Display according to the size of the window.



It becomes effective when the setting of calibration is completed. Turn grids 1.0cm apart on and off.



Return to image display from all cycle display.





Set the image display on / off. The Correct Path is automatically turned off when the image display is turned on.



During movie playback, the upper jaw is drawn at a fixed marker position. The image display is automatically turned off when the Correct Path is turned on.



.....View All Path

Display all marker positions of each tracked frame. When turned on, the image display is automatically turned off and Correct Path is turned on. It will return to the image display on Stand-By.

torrection

on Correction

When display all cycles, display the marker position with the upper jaw fixed.

		Comment	2003/7/29 🗸	Shooting Date	Frontal \sim	Direction
--	--	---------	-------------	---------------	----------------	-----------

Direction Select the shooting direction.

Shooting Date ... Designation of shooting date

CommentEnter any comment

2-2. Select Cycle / Reference



Displays one cycle of the trajectory connecting the markers and the cycle of the average.



Selected Cycle...... The selection cycle is enclosed in yellow on the left screen. The currently selected cycle is displayed in yellow on the screen on the right.

Effective Cycle...... What is checked on the left screen is the valid cycle. The effective cycle is reflected in the average cycle. The cycle being checked is displayed in green on the screen on the right.

Average Cycle...... Average value of valid cycles checked. This average cycle is reflected in "Compare Data". It is displayed in red on the screen on the right.

- **SD Mode.....**Display in standard deviation mode.
- All Valid..... Activate all cycles.
- All Invalid..... Disable all cycles.

Speed Cycle.....Display the speed cycle dialog.

Print.....You can print the displayed data. (Paper by default is set to the side of A4) **CSV Output.....**Outputs the numerical data of the table as a CSV file. This data can be edited with spreadsheet software.

* The length can be measured by dragging the mouse on the center screen.



Speed cycle dialog. It is a figure that represents the movement of the jaw in speed.

The content of each button is the same as [Cycle List].

Cycle List..... Display Cycle List dialog.

2-3. Compare Path

It is used when comparing with past analysis data.



Display "Path Comparison Screen" dialog.



Select [2 Windows] or [4 Windows], select the screen you want to insert (the screen indicating that it is selected is highlighted in yellow), and click [Insert Data] or double-click the paste area. Select comparison data dialog is displayed.

s	Select Comparison Data X											
	* Select Data to be Compared from List and Press [OK] button.											
	Analysis	Date	Age	Comment								
	Frontal	2019/ 7/30	39y2m									
	Frontal	2003/ 7/29	23y2m									
	Sagittal	2003/ 7/29	23y2m									
	Horizontal	2003/ 7/29	23y2m									
	<			>								
				OK Cancel								

Select the data to compare and click [OK].

Analysis data (average cycle) is displayed on the selected screen.



* Data not set for tracking and calibration is not displayed.

* Click [Clear] to clear the selected data.

Print.....You can print the displayed data. (Paper by default is set to the side of A4)

3. 3D Viewer

Explanation of 3D Viewer



(1)Tool bar

Exit..... Exit 3D viewer.

Save......The saved data is displayed in the analysis list.

Play..... Play the marker movement with the synthesized data.

Average...Show / Hide of average cycle.

Actual.....Show / Hide of actual cycle.

Scale.....Show / Hide of numerical scale.

Settings.... Set the background color etc.

(2) 3D display unit

Jaw movement is displayed in 3D, and you can change the viewpoint or direction by dragging the mouse.

4. Tool bar and Menu bar

4-1. Tool bar

SettingSetting

Click [Setting] to display the Preferences dialog.

Preferences	×
Set Cycle Detection • For High-Rate O For Low-Rate	Reset to Initial Value
Set Marker Detection Determination Rate 80 🔹 % Sensitivity (Color Tolerance) 60 💌	Reset to Initial Value
Set Display Color Open Close Path	
Average Path	
Character	Reset to Initial Value
Show Graph Scale <mark>Fit ∽</mark> Grid ◯ 1 mm ● 5 mm	Reset to Initial Value
Set Use Movie DirectX avi,mpg,mpeg	
QuickTime mov;mp4 * Enter Only the Extension Without a Period, and Separate Multiple With Ex: MPG;AVI	Reset to Initial Value Semicolons.
OK Cancel	

Set Cycle Detection ... Specify according to the video frame rate. (High for more than 15 FPS, Low for less than 15 FPS)

Set Marker Detection

 Determination Rate...Determination rate of marker detection

 If the software does not recognize the marker, lower rate.

 Sensitivity......

 Sensitivity (color tolerance)

 If the software does not recognize the marker, raise it.

Set Display Color...You can change the display color of each item in "Cycle List" "Speed Cycle" "Path Comparison Screen". Open is in the process of opening the mouth, Close is in the process of closing the mouth.

Show Graph...... Select the grid unit to display in the graph.

Set Use Movie...... Specify the type of movie shot with your digital camera etc.



.....Cycle Graph

Display the marker position corrected for the upper jaw with the tracked data.

<u>ک</u>	ycle Gra	aph									-		×	
	Print	Path P	oint											[Print] Print the graph.
35	-30 -:	25 -20	-15	-10) .	5)	5 '	10 1	52	02	53	03	
0														[Path Point]
0							k_							Show / Hide points (marker points)
							\mathbb{N}							
15														
20														
25														
30							Ŵ							
35							ÿ							
40														
	9	?												

.....Version information

Version

Display the Version Information dialog.

4-2. Menu bar

a. File

Preferences......Display the Preferences dialog.

Exit Analysis......Exit DigiGnatho application.

b. View

Tool Bar.....Set Show / Hide of the tool bar.

Status Bar.....Set Show / Hide of the status bar.

Display side by side at the bottom of the window.

c. Help

About DigiGnatho.... Displays the "About DigiGnatho" dialog.

* Attention:

- 1. Unauthorized reproduction of all or part of the contents of this manual is prohibited.
- 2. The contents of this program and this document are subject to change without notice for functional improvement.
- 3. The contents of this document have been prepared with every effort. Please contact us if you notice any suspicious points, errors, or omissions.
- 4. DirectX and Windows are registered trademarks of Microsoft Corporation.
- 5. Company names and product names are trademarks or registered trademarks of their respective owners.
- 6. Please note that we are not responsible for the consequences of using this program.

Developer

 Rise Corporation

 〒983-0852

 3-9-15, Tsutsujigaoka, Miyagino-ku, Sendai, Miyagi, Japan

 Tel :+81-22-295-7321

 Fax : +81-22-295-7328

 URL : http://www.risecorp.co.jp/en/